



# TSC22D1 Antibody, HRP conjugated

|                            |                                                                                                                                                                                                                     |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA623011LB01HU                                                                                                                                                                                                  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                                                                                       |
| <b>Uniprot No.</b>         | Q15714                                                                                                                                                                                                              |
| <b>Immunogen</b>           | Recombinant Human TSC22 domain family protein 1 protein (1-144AA)                                                                                                                                                   |
| <b>Raised In</b>           | Rabbit                                                                                                                                                                                                              |
| <b>Species Reactivity</b>  | Human                                                                                                                                                                                                               |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                                               |
| <b>Form</b>                | Liquid                                                                                                                                                                                                              |
| <b>Conjugate</b>           | HRP                                                                                                                                                                                                                 |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4                                                                                                                                    |
| <b>Purification Method</b> | >95%, Protein G purified                                                                                                                                                                                            |
| <b>Isotype</b>             | IgG                                                                                                                                                                                                                 |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                          |
| <b>Alias</b>               | TSC22 domain family protein 1 (Cerebral protein 2) (Regulatory protein TSC-22) (TGFB-stimulated clone 22 homolog) (Transforming growth factor beta-1-induced transcript 4 protein), TSC22D1, KIAA1994 TGFB1I4 TSC22 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                                                                                                                |
| <b>Research Area</b>       | Epigenetics and Nuclear Signaling                                                                                                                                                                                   |
| <b>Target Names</b>        | TSC22D1                                                                                                                                                                                                             |