





TSGA13 Antibody

Product Code	CSB-PA842737LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96PP4
Immunogen	Recombinant Human Testis-specific gene 13 protein (1-275AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Testis-specific gene 13 protein, TSGA13
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TSGA13


Image



Western blot
All lanes: TSGA13 antibody at 4µg/ml + A549 whole cell lysate
Secondary
Goat polyclonal to rabbit at 1/10000 dilution
Predicted band size: 32 kDa
Observed band size: 32 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA842737LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA842737LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.