



# TSGA13 Antibody, HRP conjugated

Product Code	CSB-PA842737LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96PP4
Immunogen	Recombinant Human Testis-specific gene 13 protein (1-275AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Testis-specific gene 13 protein, TSGA13
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TSGA13