



TSPAN1 Antibody, Biotin conjugated

Product Code	CSB-PA025146LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60635
Immunogen	Recombinant Human Tetraspanin-1 protein (110-211AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetraspanin-1 (Tspan-1) (Tetraspan NET-1) (Tetraspanin TM4-C), TSPAN1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TSPAN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.