



TSPAN15 Antibody, Biotin conjugated

Product Code	CSB-PA025152LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95858
Immunogen	Recombinant Human Tetraspanin-15 protein (115-235AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetraspanin-15 (Tspan-15) (Tetraspan NET-7) (Transmembrane 4 superfamily member 15), TSPAN15, NET7 TM4SF15
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TSPAN15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.