



TSPAN16 Antibody, HRP conjugated

Product Code	CSB-PA892335LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKR8
Immunogen	Recombinant Human Tetraspanin-16 protein (116-160AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetraspanin-16 (Tspan-16) (Tetraspanin TM4-B) (Transmembrane 4 superfamily member 16), TSPAN16, TM4SF16
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TSPAN16
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.