



TSPAN5 Antibody, Biotin conjugated

Product Code	CSB-PA025163LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62079
Immunogen	Recombinant Human Tetraspanin-5 protein (114-232AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetraspanin-5 (Tspan-5) (Tetraspan NET-4) (Transmembrane 4 superfamily member 9), TSPAN5, TM4SF9
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TSPAN5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.