



TSPAN6 Antibody, HRP conjugated

Product Code	CSB-PA025164LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43657
Immunogen	Recombinant Human Tetraspanin-6 protein (115-208AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetraspanin-6 (Tspan-6) (A15 homolog) (Putative NF-kappa-B-activating protein 321) (T245 protein) (Tetraspanin TM4-D) (Transmembrane 4 superfamily member 6), TSPAN6, TM4SF6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TSPAN6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.