





EIPR1 Antibody

Product Code	CSB-PA681508LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53HC9
Immunogen	Recombinant Human EARP and GARP complex-interacting protein 1 protein (201-387AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	EARP and GARP complex-interacting protein 1 (Endosome-associated recycling protein-interacting protein) (Golgi-associated retrograde protein-interacting protein) (Tumor-suppressing STF cDNA 1 protein) (Tumor-suppressing subchromosomal transferable fragment candidate gene 1 protein), EIPR1, TSSC1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	EIPR1

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA681508LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA681508LA01HU at dilution of 1:100