





# EIPR1 Antibody

<b>Product Code</b>	CSB-PA681508LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q53HC9
<b>Immunogen</b>	Recombinant Human EARP and GARP complex-interacting protein 1 protein (201-387AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	EARP and GARP complex-interacting protein 1 (Endosome-associated recycling protein-interacting protein) (Golgi-associated retrograde protein-interacting protein) (Tumor-suppressing STF cDNA 1 protein) (Tumor-suppressing subchromosomal transferable fragment candidate gene 1 protein), EIPR1, TSSC1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	EIPR1

## Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA681508LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA681508LA01HU at dilution of 1:100

---

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.