



TTC12 Antibody, Biotin conjugated

Product Code	CSB-PA862040LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H892
Immunogen	Recombinant Human Tetratricopeptide repeat protein 12 protein (1-260AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetratricopeptide repeat protein 12 (TPR repeat protein 12), TTC12
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TTC12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.