



TTC30B Antibody, HRP conjugated

Product Code	CSB-PA818703LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N4P2
Immunogen	Recombinant Human Tetratricopeptide repeat protein 30B protein (261-600AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetratricopeptide repeat protein 30B (TPR repeat protein 30B), TTC30B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TTC30B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.