



TTC38 Antibody, Biotin conjugated

Product Code	CSB-PA719094LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5R3I4
Immunogen	Recombinant Human Tetratricopeptide repeat protein 38 protein (2-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tetratricopeptide repeat protein 38 (TPR repeat protein 38), TTC38
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TTC38