



TUBA1A Antibody, Biotin conjugated

Product Code	CSB-PA754656LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q71U36
Immunogen	Recombinant Human Tubulin alpha-1A chain protein (30-229AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tubulin alpha-1A chain (Alpha-tubulin 3) (Tubulin B-alpha-1) (Tubulin alpha-3 chain) [Cleaved into: Detyrosinated tubulin alpha-1A chain], TUBA1A, TUBA3
Immunogen Species	Homo sapiens (Human)
Research Area	Isotype/Loading Controls
Target Names	TUBA1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.