



UBA7 Antibody

Product Code	CSB-PA025421LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41226
Immunogen	Recombinant Human Ubiquitin-like modifier-activating enzyme 7 protein (1-265AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-like modifier-activating enzyme 7 (Ubiquitin-activating enzyme 7) (D8) (Ubiquitin-activating enzyme E1 homolog), UBA7, UBE1L UBE2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	UBA7