





UBE2G2 Antibody

Product Code	CSB-PA025454LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60604
Immunogen	Recombinant Human Ubiquitin-conjugating enzyme E2 G2 protein (1-165AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-conjugating enzyme E2 G2 (EC 2.3.2.23) (E2 ubiquitin-conjugating enzyme G2) (Ubiquitin carrier protein G2) (Ubiquitin-protein ligase G2), UBE2G2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	UBE2G2


Image



Western blot
 All lanes: UBE2G2 antibody at 8µg/ml + K562 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 19, 16 kDa
 Observed band size: 19 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA025454LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA025454LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.