



UBE2G2 Antibody, Biotin conjugated

Product Code	CSB-PA025454LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60604
Immunogen	Recombinant Human Ubiquitin-conjugating enzyme E2 G2 protein (1-165AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-conjugating enzyme E2 G2 (EC 2.3.2.23) (E2 ubiquitin-conjugating enzyme G2) (Ubiquitin carrier protein G2) (Ubiquitin-protein ligase G2), UBE2G2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	UBE2G2