



UBE2O Antibody, FITC conjugated

Product Code	CSB-PA866339LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9C0C9
Immunogen	Recombinant Human (E3-independent) E2 ubiquitin-conjugating enzyme protein (10-277AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	(E3-independent) E2 ubiquitin-conjugating enzyme (EC 2.3.2.24) (E2/E3 hybrid ubiquitin-protein ligase UBE2O) (Ubiquitin carrier protein O) (Ubiquitin-conjugating enzyme E2 O) (Ubiquitin-conjugating enzyme E2 of 230 kDa) (Ubiquitin-conjugating enzyme E2-230K) (Ubiquitin-protein ligase O), UBE2O, KIAA1734
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	UBE2O
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.