



UBXN6 Antibody, Biotin conjugated

Product Code	CSB-PA874859LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZV1
Immunogen	Recombinant Human UBX domain-containing protein 6 protein (1-388AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UBX domain-containing protein 6 (UBX domain-containing protein 1), UBXN6, UBXD1 UBXDC2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	UBXN6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.