



UFC1 Antibody, FITC conjugated

Product Code	CSB-PA897483LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3C8
Immunogen	Recombinant Human Ubiquitin-fold modifier-conjugating enzyme 1 protein (1-167AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-fold modifier-conjugating enzyme 1 (Ufm1-conjugating enzyme 1), UFC1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	UFC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.