



UFSP2 Antibody, FITC conjugated

Product Code	CSB-PA885731LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NUQ7
Immunogen	Recombinant Human Ufm1-specific protease 2 protein (63-343AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ufm1-specific protease 2 (UfSP2) (EC 3.4.22.-), UFSP2, C4orf20
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	UFSP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.