





UGT2B4 Antibody

Product Code	CSB-PA025597LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06133
Immunogen	Recombinant Human UDP-glucuronosyltransferase 2B4 protein (281-481AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UDP-glucuronosyltransferase 2B4 (UDPGT 2B4) (EC 2.4.1.17) (HLUG25) (Hydoxycholeic acid-specific UDPGT) (UDPGTh-1), UGT2B4, UGT2B11
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	UGT2B4

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA025597LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA025597LA01HU at dilution of 1:100