



UNC50 Antibody, Biotin conjugated

Product Code	CSB-PA687496LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53HI1
Immunogen	Recombinant Human Protein unc-50 homolog protein (15-82AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein unc-50 homolog (Periodontal ligament-specific protein 22) (PDLs22) (Protein GMH1 homolog) (hGMH1) (Uncoordinated-like protein), UNC50, UNCL
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	UNC50
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.