



USMG5 Antibody, Biotin conjugated

Product Code	CSB-PA836231LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96IX5
Immunogen	Recombinant Human Up-regulated during skeletal muscle growth protein 5 protein (2-58AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Up-regulated during skeletal muscle growth protein 5 (Diabetes-associated protein in insulin-sensitive tissues) (HCV F-transactivated protein 2), USMG5, DAPIT HCVFTP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	USMG5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.