




UTP15 Antibody

Product Code	CSB-PA840580LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TED0
Immunogen	Recombinant Human U3 small nucleolar RNA-associated protein 15 homolog protein (234-518AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	U3 small nucleolar RNA-associated protein 15 homolog, UTP15
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	UTP15

Image



Western blot
 All lanes: UTP15 antibody at 5µg/ml
 Lane 1: Hela whole cell lysate
 Lane 2: HepG2 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 59, 37, 57 kDa
 Observed band size: 59 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.