



UTP15 Antibody, Biotin conjugated

Product Code	CSB-PA840580LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TED0
Immunogen	Recombinant Human U3 small nucleolar RNA-associated protein 15 homolog protein (234-518AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	U3 small nucleolar RNA-associated protein 15 homolog, UTP15
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	UTP15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.