



UTP23 Antibody, HRP conjugated

Product Code	CSB-PA874796LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BRU9
Immunogen	Recombinant Human rRNA-processing protein UTP23 homolog protein (1-249AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	rRNA-processing protein UTP23 homolog, UTP23, C8orf53
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	UTP23
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.