




VAT1L Antibody

Product Code	CSB-PA875723LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HCJ6
Immunogen	Recombinant Human Synaptic vesicle membrane protein VAT-1 homolog-like protein (1-419AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Synaptic vesicle membrane protein VAT-1 homolog-like (EC 1.-.-.-), VAT1L, KIAA1576
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	VAT1L

Image



Western blot


All lanes: VAT1L antibody at 0.3µg/ml + U87 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 46 kDa

Observed band size: 46 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA875723LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.