



VAT1L Antibody, Biotin conjugated

Product Code	CSB-PA875723LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HCJ6
Immunogen	Recombinant Human Synaptic vesicle membrane protein VAT-1 homolog-like protein (1-419AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Synaptic vesicle membrane protein VAT-1 homolog-like (EC 1.-.-.-), VAT1L, KIAA1576
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	VAT1L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.