



VEGFA Antibody

Product Code	CSB-PA07529A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15692
Immunogen	Recombinant Human Vascular endothelial growth factor A protein (27-191AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vascular endothelial growth factor A (VEGF-A) (Vascular permeability factor) (VPF), VEGFA, VEGF
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	VEGFA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.