



VPREB3 Antibody, HRP conjugated

Product Code	CSB-PA891724LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKI3
Immunogen	Recombinant Human Pre-B lymphocyte protein 3 protein (21-122AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pre-B lymphocyte protein 3 (N27C7-2) (Protein VPreB3), VPRED3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	VPRED3