



VPS26B Antibody, HRP conjugated

Product Code	CSB-PA676824LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q4G0F5
Immunogen	Recombinant Human Vacuolar protein sorting-associated protein 26B protein (1-336AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolar protein sorting-associated protein 26B (Vesicle protein sorting 26B), VPS26B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	VPS26B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.