



VWA5A Antibody, Biotin conjugated

Product Code	CSB-PA025954LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00534
Immunogen	Recombinant Human von Willebrand factor A domain-containing protein 5A protein (1-415AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	von Willebrand factor A domain-containing protein 5A (Breast cancer suppressor candidate 1) (BCSC-1) (Loss of heterozygosity 11 chromosomal region 2 gene A protein), VWA5A, BCSC1 LOH11CR2A
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	VWA5A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.