



TMEM270 Antibody, HRP conjugated

Product Code	CSB-PA744250LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UE05
Immunogen	Recombinant Human Transmembrane protein 270 protein (206-265AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 270 (Williams-Beuren syndrome chromosomal region 28 protein), TMEM270, WBSCR28
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TMEM270
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.