



WDR61 Antibody, Biotin conjugated

Product Code	CSB-PA872416LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZS3
Immunogen	Recombinant Human WD repeat-containing protein 61 protein (2-305AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	WD repeat-containing protein 61 (Meiotic recombination REC14 protein homolog) (SKI8 homolog) (Ski8) [Cleaved into: WD repeat-containing protein 61, N-terminally processed], WDR61
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	WDR61
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.