




WDTC1 Antibody

Product Code	CSB-PA822706LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N5D0
Immunogen	Recombinant Human WD and tetratricopeptide repeats protein 1 protein (201-420AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	WD and tetratricopeptide repeats protein 1, WDTC1, KIAA1037
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction

Image



Western blot

All lanes: WDTC1 antibody at 1.5µg/ml

Lane 1: A431 whole cell lysate

Lane 2: Mouse liver tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 76, 70, 68, 69 kDa

Observed band size: 76 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.