



WDTC1 Antibody, FITC conjugated

Product Code	CSB-PA822706LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N5D0
Immunogen	Recombinant Human WD and tetratricopeptide repeats protein 1 protein (201-420AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	WD and tetratricopeptide repeats protein 1, WDTC1, KIAA1037
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction