



WDYHV1 Antibody, Biotin conjugated

Product Code	CSB-PA850306LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96HA8
Immunogen	Recombinant Human Protein N-terminal glutamine amidohydrolase protein (1-205AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein N-terminal glutamine amidohydrolase (EC 3.5.1.122) (Protein NH2-terminal glutamine deamidase) (N-terminal Gln amidase) (Nt(Q)-amidase) (WDYHV motif-containing protein 1), WDYHV1, C8orf32 NTAQ1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	WDYHV1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.