



Wnt3 Antibody, Biotin conjugated

Product Code	CSB-PA026135LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17553
Immunogen	Recombinant Mouse Proto-oncogene Wnt-3 protein (22-355AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proto-oncogene Wnt-3 (Proto-oncogene Int-4), Wnt3, Int-4 Wnt-3
Immunogen Species	Mus musculus (Mouse)
Research Area	Neuroscience
Target Names	Wnt3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.