





XAGE1A Antibody

Product Code	CSB-PA872552LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HD64
Immunogen	Recombinant Human X antigen family member 1 protein (1-81AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	X antigen family member 1 (XAGE-1) (Cancer/testis antigen 12.1) (CT12.1) (G antigen family D member 2), XAGE1A; XAGE1B; XAGE1C; XAGE1D; XAGE1E, GAGED2 XAGE1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	XAGE1A

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA872552LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA872552LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.