





YIPF4 Antibody

Product Code	CSB-PA861129LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BSR8
Immunogen	Recombinant Human Protein YIPF4 protein (1-113AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:1000-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein YIPF4 (YIP1 family member 4), YIPF4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	YIPF4

Image



Western blot
 All lanes: YIPF4 antibody at 2.5µg/ml
 Lane 1: Jurkat whole cell lysate
 Lane 2: HepG2 whole cell lysate
 Lane 3: Hela whole cell lysate
 Lane 4: A431 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 28 kDa
 Observed band size: 28 kDa



Immunofluorescent analysis of HepG2 cells using CSB-PA861129LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.