



YIPF5 Antibody, Biotin conjugated

Product Code	CSB-PA842626LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969M3
Immunogen	Recombinant Human Protein YIPF5 protein (1-124AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein YIPF5 (Five-pass transmembrane protein localizing in the Golgi apparatus and the endoplasmic reticulum 5) (Smooth muscle cell-associated protein 5) (SMAP-5) (YIP1 family member 5) (YPT-interacting protein 1 A), YIPF5, FINGER5 YIP1A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	YIPF5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.