



YIPF6 Antibody, Biotin conjugated

Product Code	CSB-PA836203LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96EC8
Immunogen	Recombinant Human Protein YIPF6 protein (2-84AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein YIPF6 (YIP1 family member 6), YIPF6
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	YIPF6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.