



ZBED6CL Antibody, FITC conjugated

Product Code	CSB-PA842669LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96FA7
Immunogen	Recombinant Human ZBED6 C-terminal-like protein (1-236AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ZBED6 C-terminal-like protein, ZBED6CL, C7orf29
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ZBED6CL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.