



ZBTB18 Antibody, FITC conjugated

Product Code	CSB-PA026604LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99592
Immunogen	Recombinant Human Zinc finger and BTB domain-containing protein 18 protein (62-364AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger and BTB domain-containing protein 18 (58 kDa repressor protein) (Transcriptional repressor RP58) (Translin-associated zinc finger protein 1) (TAZ-1) (Zinc finger protein 238) (Zinc finger protein C2H2-171), ZBTB18, RP58 TAZ1 ZNF238
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZBTB18
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.