



ZBTB43 Antibody, Biotin conjugated

Product Code	CSB-PA026335LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43298
Immunogen	Recombinant Human Zinc finger and BTB domain-containing protein 43 protein (165-467AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger and BTB domain-containing protein 43 (Zinc finger and BTB domain-containing protein 22B) (Zinc finger protein 297B) (ZnF-x), ZBTB43, KIAA0414 ZBTB22B ZNF297B
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZBTB43
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.