



ZC3H10 Antibody

Product Code	CSB-PA846649LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96K80
Immunogen	Recombinant Human Zinc finger CCCH domain-containing protein 10 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger CCCH domain-containing protein 10, ZC3H10, ZC3HDC10
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ZC3H10