



ZC3H12A Antibody, HRP conjugated

Product Code	CSB-PA701732LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5D1E8
Immunogen	Recombinant Human Endoribonuclease ZC3H12A protein (1-260AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endoribonuclease ZC3H12A (EC 3.1.-.-) (Monocyte chemotactic protein-induced protein 1) (MCP-induced protein 1) (MCPIP-1) (Regnase-1) (Reg1) (Zinc finger CCCH domain-containing protein 12A), ZC3H12A, MCPIP MCPIP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ZC3H12A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.