



ZC3H7A Antibody, Biotin conjugated

| | |
|---------------------|---|
| Product Code | CSB-PA818245LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8IWR0 |
| Immunogen | Recombinant Human Zinc finger CCCH domain-containing protein 7A protein (1-167AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Zinc finger CCCH domain-containing protein 7A, ZC3H7A, ZC3H7 ZC3HDC7 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | ZC3H7A |