



ZCCHC17 Antibody, HRP conjugated

Product Code	CSB-PA878834LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NP64
Immunogen	Recombinant Human Nucleolar protein of 40 kDa protein (1-241AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleolar protein of 40 kDa (pNO40) (Pnn-interacting nucleolar protein) (Putative S1 RNA-binding domain protein) (PS1D protein) (Zinc finger CCHC domain-containing protein 17), ZCCHC17, PS1D
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZCCHC17
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.